



Pernštýnská 116,533 41 Lázně
Bohdaneč, Czech Republic
tel&fax.: +42-40-6921885
sale@rmi.cz, www.rmi.cz

Hard Coal Testing



Dilatometer RMI AA04/08

ISO 349:1975

ISO 8264: 1989

- 4 dilatometric tubes-retorts
- Very easy operation
- In-built Intel Pentium industry standard PC computer
- Passive Cooling of Electronic
- Robust Construction
- Fast Cooling of Furnace



Hard Coal Testing




Dilatation - Principle

A test piece, in the form of a pencil, prepared from powdered coal is heated at a constant rate 3K/min in a steel retort positioned in a furnace.

The change in level of a piston resting upon the test piece is observed continuously, and a record is produced which is characteristic of the swelling properties of the coal.

Hard Coal Testing

DILATATION - The measurement of the swelling of hard coal using a dilatometer



Hard Coal Testing



Dilatation - Definitions

softening temperature; temperature of initial contraction

the temperature at which the downward movement of the dilatometer piston is 0,5 mm

temperature of maximum contraction

the temperature at which the dilatometer piston reaches its lowest point

resolidification temperature; temperature of maximum dilatation

the temperature at which the dilatometer piston reaches its highest point

maximum contraction

the maximum downward movement of the dilatometer piston, measured from the zero point and expressed as a percentage of the initial test piece length

maximum dilatation

the maximum upward movement of the dilatometer piston after contraction, measured from the zero point and expressed as a percentage of the initial test piece length

Hard Coal Testing



Dilatation - Definitions

repeatability

the maximum acceptable difference between two determinations which are carried out in the same laboratory, by the same operator with the same apparatus, on test pieces prepared from the same test sample and tested simultaneously in two different retorts during the same heating cycle or separately in the same retort during different heating cycles

reproducibility

the maximum acceptable difference between the means of two determinations which are carried out in each of two laboratories, on representative portions taken from the same gross sample, after the last stage of sample preparation

Hard Coal Testing



Dilatation – Preparation of Sample

Sampling – min. weight – see ISO
Size reduction – 0.16 mm

- **New ISO** - passing 212 μm test sieve 100 %;
passing 125 μm test sieve 70 % to 60 %;
pass 63 μm test sieve 40 % to 30 %.

Storage of sample – oxidation problem

If necessary, store the test sample in an inert atmosphere in a sealed glass phial

Hard Coal Testing

Dilatation – Preparation of Sample

Preparation of test piece

10 g of sample + app. 1 ml water – mixing

- mould, mallet (hammer) and press

- test piece with dimensions

 - 60 mm $\pm$ 0,25 mm

- 60 mm $\pm$ 0,5 mm (new ISO)

Hard Coal Testing



Dilatation – Procedure

1. Place a test piece narrow end uppermost in the dilatometer retort and insert the piston into the retort so that it rests on the test piece.
2. Stabilize the furnace temperature at 30 °C below the expected softening temperature or, if this is unknown, at the appropriate temperature given in Table 2 or, if the volatile matter content is unknown, at 265 °C.
3. Insert the retort containing the test piece and piston in the appropriate hole in the furnace block and allow 10 min for the system to regain the stabilized charging temperature. Attach the recording mechanism to the piston during this period and start the test.

Hard Coal Testing

Dilatation – Procedure

Table 2 — Furnace temperature

Volatile matter content of coal (dry, ash-free basis) % (m/m)	Stabilized charging temperature °C
38,1 and above	265
28,1 to 38,0	295
18,1 to 28,0	325
18,0 and below	355

Hard Coal Testing

Dilatometr Furnace

A furnace for heating of four retorts. Starting temperature settable, end temperature 550 °C at a rate of 3 K/min. The furnace comply all specifications of all ISO norms. Low mass special design of furnace with excellent temperature uniformity and fast active cooling.

Temperature controller – PID auto adaptive regulation. Max. difference from linearity 1°C.



Hard Coal Testing

Sensors

- unique design of sensors
- LVDT types
- robust construction
- excellent stability
- easy maintenance



Specifications:

Standardised Sample Height

60 mm - see ISO 349

Measuring Range

200 mm


Sensitivity of the LVDT

< 0,5 mm


Hard Coal Testing



Electronics



Hard Coal Testing



Hard Coal Testing



Software and connections

- DR DOS
- all software is in permanent memory – can be downloaded – upgraded
- excellent stability
- very easy to use
- designed for „routine low level labs“
- data storage – standard Compact Flash Cards
- data export – via CF cards
- data – Excell compatible
- COM printer port
- Ethernet – as option

Hard Coal Testing



AA 04 – one independent unit


AA 04/08 – two independent units with better price

Advantages

- much more better flexibility
- speed of test – large furnaces needs long time for cooling !! (two or more times)

We are FAST and FLEXIBLE


Hard Coal Testing



Hard Coal Testing

Gieseler Plastometer PLAST PM04

ASTM D2639-90




Hard Coal Testing



Gieseler Plastometer PLAST PM04

ASTM D2639-90

- One retort
- Very easy operation
- In-built Intel Pentium industry standard PC computer
- Passive Cooling of Electronic
- Robust Construction
- Fast Cooling of Furnace




Hard Coal Testing



Plasticity - Principle

Measurements of the plastic properties of coals are made by applying a constant torque to a stirrer placed in a crucible into which the coal is charged. The crucible is immersed in a bath (PbSn solder) and the temperature increased uniformly (3 K/s). The rotation of the stirrer is recorded in relation to increase in temperature.

Hard Coal Testing



Hard Coal Testing



Plasticity - Definitions

Initial softening temperature

The temperature at which dial movement or electronic readout indicates a stirring shaft movement of 1.0 dial divisions per minute (dd/min), with continued indication of movement of at least 1.0 dd/min thereafter.

Maximum fluidity temperature

The temperature at which stirring shaft rotation reached the maximum rate.

Plastic range

The difference between the initial softening temperature and the solidification temperature.

Final fluidity temperature (new ISO)

The temperature at which stirrer rotation rate last reached 1 dd/min.

Solidification temperature

The temperature at which the stirring shaft stops.

Maximum fluidity

The maximum rate of rotation for the stirring shaft in dial division per minute.

Hard Coal Testing



Plasticity - Definitions

repeatability

For duplicate tests carried out at different times in the same laboratory, by the same operator with the same apparatus, on representative portions taken from the same sample after the last stage of the reduction process, all characteristic temperature points should agree within 7°C.

Likewise, maximum rates of dial drum movement of duplicate tests should agree with the repeatability values given above:

- **three results is less than 24 percent of the average**
- **the four results is less than 26 percent of the average**
- **the five results is less than 28 percent of the range, the average of all five shall be reported.**

reproducibility

Due to the inherent variance of test equipment, the values shown for Reproducibility are indicative only.

Experience with national Gieseler standards indicates that it is very difficult to assign tight limits.

Hard Coal Testing



Plasticity – Sample Preparation

- Gross sample of coal in accordance with ISO 1988 and ISO 13909
- The drying temperature should not exceed 40°C – oxidation problem
- Grinding – size reduction –see ASTM or new ISO

Hard Coal Testing



Plasticity – Sample Preparation

- Gross sample of coal in accordance with ISO 1988 and ISO 13909
- The drying temperature should not exceed 40°C – oxidation problem
- Grinding – size reduction –see ASTM or new ISO

Hard Coal Testing


Plasticity – Test Retorte

Retort crucible

Cylindrical 21.4 ± 0.1 mm inside diameter and 35 ± 0.3 mm in depth with exterior threads for joining the crucible to the barrel. The crucible shall have a 2.38 ± 0.02 mm diameter notch with an included angle of 70° in the centre of its inside base to serve as a seat for the stirrer.

Retort crucible cover

With interior threads for joining the crucible cover to the crucible and exterior threads for joining the crucible cover to the barrel. The inside diameter of the hole which accommodates the stirrer is 9.5 ± 0.1 mm



DIMENSIONS IN MILLIMETRES

FIGURE 1 — RETORT ASSEMBLY

Key

- a crucible
- b crucible cover
- c guide sleeve
- d gas exit tube
- e barrel
- f rabble arm stirrer

A diameter to fit head coupling
 B length sufficient for 75 mm immersion
 C 70° included angle
 D 1.6 mm clearance lower arm to crucible

12.7
 12.7
 3.2 ± 0.5
 3.2 ± 0.5
 35 ± 0.3
 21.4 ± 0.1
 35 ± 0.5

DIMENSIONS IN MILLIMETRES

FIGURE 1 — RETORT ASSEMBLY

Hard Coal Testing



Plasticity – Procedure

- Adjust the hysteresis brake to 101.6 g.cm $\pm$ 5 g.cm
- Mount the crucible, with the stirrer in place, in the loading device shown in Figure 5 but with the 10 kg weight removed and packing tamper raised. Charge 5.0 g of the prepared sample of coal into the crucible. Rotate the stirrer carefully with fingers in order to fill the voids under the rabble arms. Lower the packing tamper weight, making sure that the total weight is transmitted to the coal charge. Compress the charge by dropping the 1 kg weight twelve times from a height of 115 mm. For coals that are difficult to pack, it is permissible to add 1 to 3 drops of toluene on the shaft of the stirrer near the surface of the coal before the weight is added. Very dry coal may benefit by the addition of several drops of water and thorough mixing prior to packing.

Hard Coal Testing



Plasticity – Procedure


- Remove the crucible from the loading device, screw on the crucible cover taking extreme care not to disturb the position of the stirrer in the coal. Screw the crucible and stirrer into the retort assembly. The stirrer should be centred in the guide sleeve. Screw the retort assembly onto the plastometer head, making sure that the top end of the stirrer fits into the slotted end of the axle.
- Lower the assembled apparatus into the furnace until the bottom of the crucible is immersed to a depth of 75 mm in the molten solder bath maintained at a temperature of 300°C. Place the thermocouple in the bath in accordance with Clause 5.3. The heating control shall be such that the bath will re-attain the initial temperature in 10 ± 2 min after immersion of the retort. Thereafter, heat the retort at a rate of $3.0 \pm 0.1^\circ\text{C}/\text{min}$ on the overall basis. When the movement of the drum or readout from the electronic sensor has reached 1.0 dial division per minute, take readings of temperature and dial movement at 1 minute intervals. Continue readings until the dial shows no further movement.

Hard Coal Testing

Plastometr Furnace

A furnace for heating of one retorts in solder bath (Pb/Sn solder). Starting temperature settable, end temperature 550 °C at a rate of 3 ± 0.1 K/min. The furnace comply all specifications of ASTM norm. Fast active cooling.

Temperature controller – PID auto adaptive regulation.



Hard Coal Testing

Plastometr Furnace

Two temperature sensors
(heating element and crucible) type K


Internal diameter of the
furnace 65 mm

Volume bath ~400 cm³

Charge of the alloy Sn50Pb 3 kg

Operating voltage 48 Vac

Input max. 800 W



Hard Coal Testing

Head

- stirrer is loaded by constant torque 101,6 g.cm.
- magnetic clutch
- automatic calibration procedure with excellent reproducibility
- excellent stability of torque
- registration – 100 divisions per cycle – optoelectric



Hard Coal Testing



Software and connections

- DR DOS
- all software is in permanent memory – can be downloaded – upgraded
- excellent stability
- very easy to use
- designed for „routine low level labs“
- data storage – standard Compact Flash Cards
- data export – via CF cards
- data – Excell compatible
- COM printer port
- Ethernet – as option

Hard Coal Testing



Advantages

- much more better precision of calibration
- better long term stability
- better serviceability